

# Benzalkonium Chloride & Glutaraldehyde Solution

## EXOTICA SD

### Aldehyde Based High level Surgface Disinfectant

Aldehyde, quaternary ammonium compound,  
And chemically bonded formaldehyde based high level  
surface disinfectant cleaner

With universal application. Safe for  
Disinfection cleaning in risk areas.

Suitable for terminal instrument disinfection.

#### Composition:

It should be

Each 100 gm contains


Benzalkonium Chloride

(eq. to Benzalkonium Chloride Solution IP) IP: 5g

1, 6 Dihydroxy, 2-5 Dioxahexane: 11.2 g

Glutaraldehyde (Eq. to Strong Glutaraldehyde Solution IP): 5 g

Alkyl Urea Derivative: 3 g.



#### Direction For Use :

- Disinfection of precleaned surface, critical and non critical instruments: 10 ml in 1 liter of water.
- Disinfection of critical areas: 20 ml in 1 litre of water. Let the surface remain wet for 10 minutes.
- Disinfection of surgical refuge: 15 ml in 1 litre of water. Dip the refuge in the solution for 40 minutes before disposal.

#### Features

- Unique blend of three potent biocide.
- Fast action .
- Long lasting residual effect.
- Penetrates even the difficult to reach areas.
- Synergic combination of Glutaraldehyde and chemically bonded formaldehyde in water creates a potent biocide having excellent microbicidal activity.
- Only 30% of the formaldehyde is used at a time rest is released only when there is an addition of new contaminant.

#### Precaution

Do not mix with detergents or chemicals.  
Keep out of the reach of children.  
Avoid direct contact with eyes and ears.

#### Storage Condition :

**Store in dry place. Protect from light and heat.**

FOR EXTERNAL USE ONLY